
Red Light Mask Laser

User Manual

Optimal Results: Frequency and Duration of Use

For the best results Red Light Mask, follow these guidelines (do not use for more than 20 minutes per day):

M1. Glow Boost Therapy

- *Recommended Use up to:* 10 minutes at high brightness
 - *Frequency:* 4 times per week
-

M2. Clear Skin Therapy

- *Recommended Use up to:* 5- 10 minutes at high brightness
 - *Frequency:* 4 times per week
-

Warranty & Guarantee Terms

offers a one-year guarantee covering material and craftsmanship defects under normal use, with exclusions for misuse, neglect, or unauthorized repairs.

Specification Details

Product Name	Red light Facial Mask Laser
Model No.	TLM300-L
Product Size	640mm * 220mm * 5.8mm
Product Weight	340g
Light Wavelength	460nm,665nm,850nm,1064nm
Output Power	5W
Chips Quantity	236pcs
Battery	1440mAh, use time 35 minutes
Function	Light Mode + Dimmable + Timer
Blue Light	460 ± 3nm Led chips 72pcs
Red Light	665±3nm Laser chips 72pcs
Near-Infrared Light	850 ± 3nm Laser chips 72pcs
Near-Infrared Light	1064±3nm Laser chips :20pcs

Setting Up Your Mask for First Use

1. Prepare Your Skin:

Ensure your face is clean and dry.

2. Fitting the Mask:

Place the mask on your face and adjust the strap for comfort.

3. Turning On:

Press and hold the power button.

4. Starting Treatment:

Press the mode button to select your desired therapy.

5. Enjoy the Session:

Relax and adhere to the recommended therapy times.

Cleaning the Mask

Wipe with a damp microfiber cloth.

Do not use abrasive materials or strong cleaning agents.

Safety Certifications

Our device complies with CE, RoHS, FCC, GB, IE meeting international safety standards


Charging the Mask

Connect to power using the USB-C cable.

Charging Indicators:

Red light: Charging

Green light: Fully Charged



Type C Charging port

Additional Notes.

Charging time: 3 hours


Battery life: Enjoy up to 40 minutes of use, allowing for 4 ten-minute sessions on low brightness or 3 ten-minute sessions on high brightness, per charge.

Cannot be used while charging.

When the battery is low, the indicator light on the top of the mask will turn red. Please charge it in time to avoid affecting the next use of the product.

USB Charging Instructions for Controller

1. Locate the USB charging port on the controller.
2. Connect the provided USB cable to the controller and a power source (e.g., wall adapter or computer).
3. The indicator light will flash during charging.
4. The indicator will remain solid when charging is complete.




Disconnect the USB cable after charging.

Light Therapy Modes and Customization

Adjust the timer to fit your skin's needs.
When first using the mask, use T1 for 5 minutes as a test

Recommended Timer Settings:

- T1:** 5-15 minutes
- T2:** 5-10 minutes
- T3:** 10-15 minutes



Brightness Levels:

(For each mode, use Brightness on HIGH)

- Level 1:** Low
- Level 2:** Medium
- Level 3:** High



Modes:

M1. Glow Boost Therapy (Red + Near Infrared Light)

- **Benefits:** Reduces wrinkles, fine lines, and boosts circulation.
- **Recommended Use up to:** 10 minutes, 4 times per week on high brightness.

M2. Clear Skin Therapy (Blue + Near Infrared Light)

- **Benefits:** Helps control acne, reduce inflammation, and redness.
- **Recommended Use up to:** 5- 10 minutes, 4 times per week on high brightness.
- **Recommended Use up to:**
 - 10-minute sessions: Up to 4 times per week on high brightness.
 - 15-minute sessions: Limit to 3 times per week on high brightness.

Note: For 15-minute sessions, start with three 10-minute sessions the previous week to condition your skin.


Adjustable Brightness Settings for Sensitive Skin:

Our mask offers lower brightness settings for individuals whose skin does not take well to the higher power of the mask

How to Use the Mask

Turning On/Off:

Press and hold the button for 2 seconds.




Change Modes:

Press the **M** button on the remote control to cycle through:

M1. Glow Boost Therapy
(Red+Double Near Infrared Light)

M2. Clear Skin Therapy
(Blue + Near Infrared Light)



Safety Instructions

For your safety and to ensure proper use, please read the user manual thoroughly before operating this device.



Caution:

Improper use may result in minor injury or damage to the device.



Notice:

Ensure all setup and operational procedures are followed as detailed in this manual to avoid potential harm or device damage.



Warning:

Avoid exposure to water or heat sources to prevent the risk of injury and device malfunction.



Specific Dangers:

Not suitable for users of pacemakers or similar devices.

Keep away from water to prevent malfunction.

Avoid use near heat sources like stoves or air conditioners.

Do not use if damaged to avoid risk of electric shock.

Keep out of reach of children.

Do not use the mask for more than 20 minutes daily.

encouraging
wellness naturally _____